DESMOCLEAN®

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DESMOCLEAN®

Instructions for Use

Rotating instrument for gentle curettage and planing of the root surface.

The slim, elliptical, non-cutting hexagonal head reaches all subgingival surfaces, even in hard-to-access areas, such as very deep pockets and furcations.

Features of the Desmoclean instrument:

The instrument fits into the low speed contra-angle. The maximal speed should be between 7,000 and 10,000 rpm.

The slim, elliptical shape reaches even into narrow approximal spaces and deep pockets.

The hexagonal head has no cutting edges, so the surface is cleaned, but not damaged. Subgingival concrements can also be removed easily.

Visual evidence by SEM


Comparative studies show the root surface treated with Desmoclean, scaler, diamond finishing bur 40 µm, diamond finishing bur 15 µm.

Desmoclean (50:1)

Scaler (51:1)

Diamond finishing bur 40 µm (51:1)

Diamond finishing bur 15 µm (50:1)


Important:

Sterilize before first use. Operate Desmoclean with only light pressure and water spray at 7,000 – 10,000 rpm. One instrument should be used for 30 – max. 50 patients, depending on the intensity of the treatments.

Summary

Desmoclean is ideal for

- curettage, open and closed
- concrement removal
- planing of tooth surfaces after removal of brackets or splinting

Shape I

1:1 = 23 mm

Shape II

1:1 = 29 mm

Shape III

1:1 = 26 mm

Shape IV

1:1 = 29 mm

Orders and Deliveries

Introductory Package:

6 instruments per shape = 24 instruments REF 503 290

Single Shape Package: 12 instruments per shape

Desmoclean I	= 12 instruments	REF 503 291
Desmoclean II	= 12 instruments	REF 503 292
Desmoclean III	= 12 instruments	REF 503 293
Desmoclean IV	= 12 instruments	REF 503 294

Combination Package: 6 each instruments of 2 shapes

Desmoclean I + II	= 12 instruments	REF 503 295
Desmoclean III + IV	= 12 instruments	REF 503 296
Desmoclean I + III	= 12 instruments	REF 503 297
Desmoclean II + IV	= 12 instruments	REF 503 298

Valid as of: 12/2012

Two head sizes

Two shaft lengths

= 4 instruments

Manufacturer Information on Reprocessing of Medical Devices according to DIN EN ISO 17664:2004

Manufacturer: HAGER & WERKEN GmbH & Co. KG
Ackerstraße 1
47269 Duisburg, Germany
T +49 (203) 99269-0
F +49 (203) 299283
www.hagerwerken.de

Product: Desmoclean
REF 503 265 - REF 503 268
REF 503 291 - REF 503 298

Warning notices: no special notices

Limitation of reprocessing: Frequent reprocessing has only got a small impact on this instrument. The instrument shelf-life is mainly determined by wear and damage by use.

Place of Handling: Clean surface with a disposable cloth or paper towel.

Storage and Transport: No special requirements.

Preparations for Cleaning: None

Cleaning: Only use suitable cleansers, carefully following the instructions of the respective manufacturer. After cleaning, remove cleanser thoroughly under running water. Do not exceed 93 °C in mechanical cleaning.

Disinfection: Only use suitable disinfecting solutions, carefully following the instructions of the respective manufacturer. Remove disinfecting solution thoroughly under running water. Do not exceed 93 °C in mechanical disinfection.

Packaging: Standardized packing material for sterilisation can be used (recommendation).

Sterilization: Vapor sterilization at 134 °C for 5 minutes or vapor sterilization at 121 °C for 20 minutes.

Control/Functional Check: Sight check on damages, wear, deformation.

Storage: No special requirements.

The a.m. instructions have been validated as SUITABLE for preparation of a medical device and its re-use by the medical device manufacturer. It is the reprocessor's responsibility that the actual reprocessing achieves the required results with the equipment and materials applied and personnel involved in the reprocessing site. Normally, validation and routine control are necessary. Furthermore, the reprocessor should evaluate any deviation from the provided instructions regarding efficiency and possible adverse consequences.