
900 3242

SINGLE COAT
BOND**Universal light-curing dental adhesive****Adhesivo dental universal fotocurable****Adhésif dentaire universel photopolymérisable****Universeller, lichthärtender Haftvermittler****Adesivo dentale universale fotopolimerizzabile****Sistema adesivo fotopolimerizável**

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EN

Definition

CyberTech Single Coat Bond is a light-curing, multifunctional single-component bonding agent for the adhesive restoration technique.

Composition

Methacrylates
Photoinitiators
Methacrylated polyalkenoate

Indications

CyberTech Single Coat Bond is indicated for:

- Direct filling using composite and compomer
- Bonding of indirect restorations (inlays, onlays, veneers and crowns) made of ceramic or composite
- Bonding of composite to ceramic, composite, metal and amalgam
- Dentin sealing

Contra-indications

Proven hypersensitivities to ingredients of CyberTech Single Coat Bond.

Side effects
The ingredients of CyberTech Single Coat Bond can cause sensitivity in predisposed persons.

Interactions

Avoid the use of phenolic containing materials and other substances which inhibit polymerization (e.g. zinc oxide/eugenol) in contact with CyberTech Single Coat Bond.

Applications

- Direct filling
- Bonding of inlays, onlays, veneers and crowns
- Bonding to other materials

1. Cavity preparation

According to the principle of the adhesive restoration technique.

For pulp protection, a small amount of hard curing calcium hydroxide cements should be placed close to the pulp. Finish the enamel again with a diamond bur. For isolation, the use of a rubber dam is recommended.

2. Conditioning

Apply phosphoric acid gel on enamel and dentin according to the manufacturer's instructions for use. Afterwards rinse thoroughly with water ensuring that all residues have been removed. Dry the cavity with oil-free air.

Important: Do not desiccate dentin.

3. Application

Dispense one drop CyberTech Single Coat Bond into a mixing well and apply it to the cavity using a disposable dentalbrush. Rub the surface for 20 s and gently dry with oil-free air. Light cure for 30 s (halogen or LED curing light with a light intensity ≥ 800mW/cm²). At lower light intensity, exposure time is extended accordingly.

Note: Avoid pooling of adhesives. This can impair the accuracy of fit.

4. Composite, compomer and other materials

Procedure according to the manufacturer's instructions for use.

Important: For the use of CyberTech Single Coat Bond must be ensured, that sufficient incident light penetrates through the complete restoration. According to the manufacturer's specifications.

Emergency measures, symptoms, antidotes

In case of direct contact with the oral mucosa, simply rinse with water. In case of contact with eyes, rinse thoroughly with water. Get medical attention for assessment and further treatment.

Notes

Only supplied to dentists and dental technicians or upon their instructions.

Keep out of the reach of children!

Only dispose of completely empty packages.

Shelf life and labelling

Expiry date and [LOT] number can be found on the containers and outer packaging.

Do not use the product beyond its expiry date.

Storage

Store at 4–23°C / 39–73°F.

Do not expose to direct sunlight or other heat sources.

Hygiene tips

For optimal hygienic conditions of using CyberTech Single Coat Bond the following instructions are recommended:

- Disposable dentalbrushes are for single patient use only.
- Close the bottles immediately after use. Clean with standard disinfectants.

Caution

Federal law restricts this device to sale by or on the order of a dentist.

Date of issue of instructions for use

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