TRAY AND INSTRUMENT CLEANER CONCENTRATE

Fluid concentrate for cleaning as part of the correct preparation of medical products

- **Powerful cleaning**
- **Particularly gentle on materials**
- **Ideal for use in ultrasonic bath**

PRODUCT DESCRIPTION:

CyberTech Tray and instrument cleaner concentrate is a concentrate for cleaning of impression trays, mixing aids and dental instruments. Because of its mild, alkaline composition, the innovative formulation offers very good cleaning performance while being gentle on materials. CyberTech Tray and instrument cleaner concentrate is suitable for cleaning all metals resistant to corrosion and sensitive materials such as anodised aluminium. CyberTech Tray and instrument cleaner concentrate is also suitable for use in ultrasonic baths.

AREAS OF APPLICATION:

Liquid concentrate for cleaning impression trays, mixing accessories and dental instruments.

CONCENTRATION:

Please dilute the concentration to the required solution by adding water.

	5 %	6 %	7 %	8 %	9 %	10 %
1 L*	50 ml	60 ml	70 ml	80 ml	90 ml	100 ml
2 L*	100 ml	120 ml	140 ml	160 ml	180 ml	200 ml
4 L*	200 ml	240 ml	280 ml	320 ml	360 ml	400 ml
8 L*	400 ml	480 ml	560 ml	640 ml	720 ml	800 ml
10 L*	500 ml	600 ml	700 ml	800 ml	900 ml	1000 ml

APPLICATION:

Use in soaking bath: Roughly pre-wash soiled impression trays and instruments and, depending on level of soiling, place in a 5 - 10% application solution (50 to 100ml cleaning concentrate to 1 litre of application solution). Leaving the solution to work overnight provides simple handling and optimal cleaning performance. Then carefully remove any remaining or not fully-dissolved dirt particles with water and a suitable brush from the impression trays and instruments. After thorough rinsing, continue to prepare the instruments.

Use in ultrasonic bath: Depending on level of soiling, treat the impression trays and instruments for 15 – 30 mins in a 10% application solution (100 ml cleaning concentrate to 1 litre of application solution) in an ultrasonic bath. For particularly stains, repeat the process. Then thoroughly rinse the cleaned items.

English

COMPOSITION:

Anionic surfactants 5 - 15%.

SPECIAL INFORMATION:

Sulfonic acid, C14-17-sec-alkan-, sodium salt; glycolic acid. Causes skin irritation. Causes serious eye damage. Wear protective gloves/protective clothing/eye protection/face protection. IF ON SKIN: Wash with plenty of water. If skin irritation occurs: Get medical advice/attention. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. Immediately call a POISON CENTER/doctor.

PRODUCT STATUS/LISTINGS:

Conform to the EU guideline 93/42/EC for medical devices.

AVAILABLE DELIVERY FORMS

Itemnumber	Packaging
9882809	1L Dosing bottle

Use disinfectants safely. Always read the label and product information before use.

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