INSTRUMENT DISINFECTION CONCENTRATE

**Economical, materialfriendly,
active-cleaning for routine
clinical use**

- **Aldehyde-free preparation with strong cleaning power**
- **High degree of material compatibility**
- **Pleasant scent**

PRODUCT DESCRIPTION:

CyberTech Instrument disinfection concentrate is an economical, material friendly, active-cleaning instrument disinfectant for routine clinical use. Thanks to the exceptional antimicrobial efficacy a high degree of personal protection is assured during instrument (re-)processing. CyberTech Instrument disinfection concentrate is characterised by a pleasant smell thanks to the aldehyde-free formulation and is suitable for use on medical instruments of all kinds.

AREAS OF APPLICATION:

Disinfection and cleaning of endoscopes and (re-) processing of all types of medical instruments.

SPECTRUM OF ACTION:

Bactericidal (incl. MRSA), yeasticidal, virus-inactivating (HBV/HIV/HCV viruses).

CONCENTRATION:

1% – 60 mins.; 2% – 30 mins.; 3% – 15 mins.

APPLICATION:

Reprocessing of all types of non-hinged instruments disinfecting cleaning by immersion bath and disinfecting by immersion bath. Suitable for use in ultrasonic baths. Disinfectant cleaning in the immersion bath procedure: After preliminary cleaning, submerge the instruments in the CyberTech Instrument disinfection solution and clean them mechanically. All surfaces and cavities must be wetted completely by the working solution. After the contact time has elapsed, rinse thoroughly with demineralised water.

Disinfectant cleaning in the ultrasonic bath procedure: CyberTech Instrument disinfection concentrate can be used to improve cleaning performance in all stainless steel ultrasonic baths. Prepare CyberTech Instrument disinfection concentrate working solution. The temperature of the working solution should not be higher than 50°C. After the contact time has elapsed, rinse the medical devices thoroughly with demineralised water.

Dry and sterilise critical instruments using a validated procedure.

APPLICATION INFORMATION:

Appropriate personal protective equipment (PPE) must be worn when handling CyberTech Instrument disinfection concentrate. For further information about occupational health and safety, see www.bgw-online.de. Before using the CyberTech Instrument disinfection concentrate, clean all instruments and baths thoroughly to remove chemical residues. Instrument baths must be cleaned and disinfected at least once every working day, preferably with a disposable wipe soaked in alcohol.

COMPOSITION:

100 g contains: 3.75 g N,N-didecyl-n-methylpoly(oxyethyl)ammoniumpropionate; 2.5 g guanidine, N,N"-1,3-propanediylbis-, N-coco alkyl derivs., diacetates.

SPECIAL INFORMATION:

N,N-didecyl-n-methyl-poly(oxyethyl)ammoniumpropionate; alcohol ethoxylate C9-C11; guanidine, N,N"-1,3-propanediylbis-, N-coco alkyl derivs., diacetates. Can cause severe skin burns and eye damage. Very toxic to aquatic life with long lasting effects. Wear protective gloves/protective clothing/eye protection/face protection. IF SWALLOWED: rinse mouth. Do NOT induce vomiting. IF ON SKIN (or hair): Take off immediately all contaminated clothing. Rinse skin with water/ shower. IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. IF exposed: Call a POISON CENTER or doctor/physician. Avoid release to the environment.

INSTRUMENT DISINFECTION CONCENTRATE

**Economical, materialfriendly,
active-cleaning for routine
clinical use**

- **Aldehyde-free preparation with strong cleaning power**
- **High degree of material compatibility**
- **Pleasant scent**

English

PRODUCT STATUS/LISTINGS:

Included in the current DGHM/VAH list of disinfectants.
Conform to the EU guideline 93/42/EC for medical devices.

AVAILABLE DELIVERY FORMS

Itemnumber	Packaging
9882813	2L Bottle
9882814	5L Canister

AVAILABLE DELIVERY FORMS

Itemnumber	Packaging
9796273	Dosing pump for 2L bottle
9796276	Tab/dosing device for 5L canister
9796278	Canister-key
9796272	Dosing pump for 5L Canister

Use disinfectants safely. Always read the label and product information before use. In case of visible soiling the CyberTech Instrument disinfection concentrate working solution must be changed immediately; otherwise, change it at least every working day.

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