

Moisture-Activated Temporary Filling and Sealling Material

SECTION I - PRODUCT / COMPANY IDENTIFICATION

Product Name: TEMPIT® TEMPORARY FILLING MATERIAL

Synonyms: Temporary filling material Part (Item) Number: 310060, 360053

Company Name: Centrix, Inc. Address: 770 River Road

City, State, Zip: Shelton CT, 06484

Emergency Telephone Number: 203-929-5582 Information Telephone Number: 800-235-5862

SECTION II - COMPOSITION / INFORMATION ON INGREDIENTS

Hazardous Components:	%	CAS No.	SYMBOL	R	EINECS No.
D.M.T.G (Glycol)	<12%	143-24-8			205-594-7
Ethyl Methacrylate Polymer	<3%				
Calcium Sulphate	<59%	10034-76-1			231-900-3
Zinc Oxide	<11%	1314-13-2			
Barium Sulphate	<17%	7727-43-7			
Silica (Amorphous)	<2%	7631-86-9			231-545-4

SECTION III - HAZARDS IDENTIFICATION

For potential health hazards: See section 11 below.

No other significant hazards to man or the environment are reasonably foreseen.

SECTION IV - FIRST AID MEASURES

Skin: Wash with soap and water.

Eyes: Flush with copious amounts of water. If irritation persists, seek medical attention.

Ingestion: N/A Inhalation: N/A

SECTION V - FIRE-FIGHTING MEASURES

Use extinguisher suitable to cause of fire and its surroundings.

Suitable extinguishing media: N/A
Special firefighting measures: N/A
Unusual fire and explosion hazards: N/A

Special protection equipment: N/A

SECTION VI - ACCIDENTAL RELEASE MEASURES

Transfer bulk into suitable containers for disposal. Wipe residue with damp cloth and scrape up into

suitable container for disposal.

Personal precautions: N/A Environmental precautions: N/A

Reclaiming methods: N/A

SECTION VII - HANDLING AND STORAGE

Handling: No handling problems are likely to be encountered provided normal safe practices of hygiene are applied.

Storage: Store in a cool dry place. Containers: N/A Incompatibilities: N/A Other precautions: N/A

SECTION VIII - EXPOSURE CONTROLS / PERSONAL PROTECTION

General measures: Wash hands thoroughly after contact with this product.

Limits: N/A

Respiratory protection: N/A

Eye protection: Goggles/eyeshields. Skin protection: Overalls, rubber gloves. Hands protection: Overalls, rubber gloves.

Other equipment: N/A

SECTION IX - PHYSICAL AND CHEMICAL PROPERTIES

Appearance: White stiff paste Odor: Odorless

Explosion limits: N/A Water solubility: Insoluble

Boiling point: N/A Evaporation rate: N/A Flash point: N/A

SECTION X - STABILITY AND REACTIVITY

Stability: Stable, non-reactive Conditions to avoid: N/A Materials to avoid: N/A

Hazardous decomposition products: N/A Hazardous polymerization products: N/A

SECTION XI - TOXICOLOGICAL INFORMATION

Ingestion: N/A
Eye contact: N/A
Skin contact: N/A
Inhalation: N/A

Prolonged exposure: N/A

SECTION XII - ECOLOGICAL INFORMATION

SECTION XIII - DISPOSAL CONSIDERATIONS

Dispose of in accordance with all federal, state, and local regulations.

SECTION XIV - TRANSPORT INFORMATION

SECTION XV - REGULATORY INFORMATION

Hazard Classification Symbol: N/A

Ingredients determining the label: The information is this MSDS does not constitute an assessment of workplace risk.

SECTION XVI - OTHER INFORMATION

Handling of product exclusively by dental specialists according to instructions for use.

CAUTION: PRODUCT FOR PROFESSIONAL USE ONLY

The above information is based on presently available data and to our best knowledge for handling the product under normal conditions. Any use of this product in any way not indicated on this document or using it together with any other process/procedure will be exclusively under the user's responsibility.



Making Dentistry Easier.[™]

CENTRIX INCORPORATED

C€ 0510

770 River Road, Shelton, CT 06484 USA (800) 235-5862 Tel: (203) 929-5582 Fax: (203) 929-6804



CENTRIX EUROPE

Hauptstrasse 64 B, 65719 Hofheim, Germany

Product covered by one or more of the following patents: US 4,963,093 and US D353673. Other U.S. and International patents pending or granted.

©Centrix, Inc. 2006 Date Revised 10/15/2010