Material Safety Data Sheet

May be used to comply with

OSHA's Hazard Communication Standard. 29 CFR 1900, 1200 Standard must be consulted for specific requirements. **U.S. Department of Labor**

Occupational Safety and Health Administration

(Non-Mandatory Form) Form Approved OMB No. 1218-0072

IDENTITY (As used on label and list)

Ti-Core Auto E (Base)

SECTION I

Manufacturer's Name

Essential Dental Systems, Inc.

Address (Number, Street, City, State, Zip)

89 Leuning Street

South Hackensack, NJ 07606

Emergency Telephone Number

Telephone Number for Information

(201) 487-9090

Date Prepared March 1, 2009

SECTION II - Hazardous Ingredients/Identity Information

Hazardous Components Other Limits

(Specific Chemical Identity) OSHA PEL ACHIH TIV Recommended %Optional

AROMATIC AND ALIPHATIC DIMETHACRYLATES

INORGANIC PIGMENTS

NN'-DIHYDROXYETHYL P-TOLUIDINE

SECTION III -Physical/Chemical Characteristics

Boiling Point Specific Gravity $(H_2O = 1)$

NE 1.

Vapor Pressure (mm Hg) Melting Point

N/A N/A

Vapor Density (AIR = 1) Evaporation Rate (Butyl Acetate = 1)

N/A N/A

Solubility in Water

INSOLUBLE

Appearance and Odor

OFF WHITE PASTE WITH CHARACTERISTIC ACRYLATE ODOR

SECTION IV - Fire and Explosion Hazard Data

Flash Point (Method Used) Flammable Limits LEL UEL N/A N/A N/A N/A

Extinguishing Media

DRY CHEMICAL FOAM, CO2, DRY CHEMICAL

Special Fire Fighting Procedures

N/A

Unusual Fire and Explosion Hazards

MATERIAL MAY BURN IF EXPOSED TO FLAME.

MSDS/Ti-Core Auto E (Base)

Page 2 of 2

SECTION V - Reactivity Data

Stability Unstable Conditions to Avoid HIGH TEMPERATURES

Stable

 \mathbf{X}

Incompatability (Materials to Avoid)

N/A

Hazardous Decomposition or Byproducts

N/A

Hazardous Polmerization May Occur Conditions to Avoid

X HIGH TEMPERATURE, STRONG OXIDIZERS AND PEROXIDES

Will not Occur

SECTION VI - Health Hazard Data

Routes of Entry: Inhalation Skin Ingestion

N/A Yes N/A

Health Hazards (Acute and Chronic)

None

Carcinogenicity: NTP IARC Monographs OSHA Regulated

NONE KNOWN

Signs and Symptoms of Exposure:

REPEATED CONTACT WITH SKIN MAY CAUSE REDDENING AND IRRITATION.

Medical Conditions Generally Aggravated by Exposure:

UNKNOWN

Emergency and First Aid Procedures:

EYES – FLUSH WITH COPIOUS AMOUNTS OF WATER. SKIN- WASH WITH SOAP AND WATER.

SECTION VII - Precautions for Safe Handling and Use:

Steps to Be Taken in Case Material is Released or Spilled:

Waste Disposal Method

DISPOSE OF IN ACCORDANCE WITH FEDERAL, STATE, AND LOCAL AUTHORITIES

Precautions to Be Taken in Handling and Storing:

STORE AT OR BELOW ROOM TEMPERATURE AVOID MOISTURE CONTACT IN

UNPOLYMERIZED STATE

Other Precautions

Protective Gloves

NONE

SECTION VIII - Control Measures

Respiratory Protection (Specify Type)

N/A

Ventilation:Local ExhaustSpecialN/AN/AN/AMechanical (General)Other

Eye Protection

RECOMMENDED RECOMMENDED Other Protective Clothing or Equipment

N/A

Work/Hygienic Practices

GENERAL SAFE PRACTICES

The data given here is based on current knowledge and experience. The data does not signify any warranty with regard to the product's properties.

Material Safety Data Sheet

May be used to comply with

OSHA's Hazard Communication Standard. 29 CFR 1900, 1200 Standard must be

consulted for specific requirements.

U.S. Department of Labor

Occupational Safety and Health Administration

(Non-Mandatory Form) Form Approved

OMB No. 1218-0072

IDENTITY (As used on label and list)

Ti-Core Auto E (Catalyst)

SECTION I

Manufacturer's Name

Essential Dental Systems, Inc.

Address (Number, Street, City, State, Zip)

89 Leuning Street

South Hackensack, NJ 07606

Emergency Telephone Number

Telephone Number for Information

(201) 487-9090

Date Prepared March 1, 2009

SECTION II - Hazardous Ingredients/Identity Information

Hazardous Components

Other Limits

(Specific Chemical Identity) OSHA PEL ACHIH TIV Recommended %Optional

AROMATIC AND ALIPHATIC DIMETHACRYLATE MONOMERS

BENZOYL PEROXIDE

2,6-DI-TER BUTYL-P-CRESOL (BHT)

SECTION III -Physical/Chemical Characteristics

Boiling Point Specific Gravity $(H_2O = 1)$

NE 1.7

Vapor Pressure (mm Hg) **Melting Point**

Vapor Density (AIR =1) **Evaporation Rate (Butyl Acetate = 1)**

N/A N/A

Solubility in Water

INSOLUBLE

Appearance and Odor

OFF WHITE PASTE WITH CHARACTERISTIC ACRYLATE ODOR

SECTION IV - Fire and Explosion Hazard Data

Flash Point (Method Used) **Flammable Limits** LEL UEL N/A N/A N/A N/A

Extinguishing Media

DRY CHEMICAL FOAM, CO₂, DRY CHEMICAL

Special Fire Fighting Procedures

N/A

Unusual Fire and Explosion Hazards

MATERIAL MAY BURN IF EXPOSED TO FLAME.

MSDS/Ti-Core Auto E (Catalyst)

Page 2 of 2

SECTION V - Reactivity Data

Stability Unstable Conditions to Avoid HIGH TEMPERATURES

Stable

X

Incompatability (Materials to Avoid)

N/A

Hazardous Decomposition or Byproducts

N/A

Hazardous Polmerization May Occur Conditions to Avoid

X HIGH TEMPERATURE, STRONG OXIDIZERS AND PEROXIDES

Will not Occur

SECTION VI - Health Hazard Data

Routes of Entry: Inhalation Skin Ingestion

N/A Yes N/A

Health Hazards (Acute and Chronic)

None

Carcinogenicity: NTP IARC Monographs OSHA Regulated

NONE KNOWN

Signs and Symptoms of Exposure:

REPEATED CONTACT WITH SKIN MAY CAUSE REDDENING AND IRRITATION.

Medical Conditions Generally Aggravated by Exposure:

UNKNOWN

Emergency and First Aid Procedures:

EYES – FLUSH WITH COPIOUS AMOUNTS OF WATER. SKIN- WASH WITH SOAP AND WATER.

SECTION VII - Precautions for Safe Handling and Use:

Steps to Be Taken in Case Material is Released or Spilled:

Waste Disposal Method

DISPOSE OF IN ACCORDANCE WITH FEDERAL, STATE, AND LOCAL AUTHORITIES

Precautions to Be Taken in Handling and Storing:

STORE AT OR BELOW ROOM TEMPERATURE AVOID MOISTURE CONTACT IN

UNPOLYMERIZED STATE

Other Precautions

NONE

SECTION VIII - Control Measures

Respiratory Protection (Specify Type)

N/A

Ventilation:Local ExhaustSpecialN/AN/AN/A

Mechanical (General) Other

Protective Gloves Eye Protection
RECOMMENDED RECOMMENDED
Other Protective Clothing or Equipment

N/A

Work/Hygienic Practices

GENERAL SAFE PRACTICES

The data given here is based on current knowledge and experience. The data does not signify any warranty with regard to the product's properties.